

Specifications

- Product** Pro-Tite™ Resin-Based Thread Sealant
- Application** Thread sealant for metals and plastics
- Indoor/Outdoor Use** Indoor or outdoor
- Storage** Cool, dry place
- Purpose** Pro-Tite™ is a resin based multi-purpose thread sealant designed to form a gap-filling seal on all types of threaded fittings and pipe. Pro-Tite™ is an easy-to use flowable paste supplied in a convenient 1.3-ounce squeeze tube, and in convenient 4 and 16-ounce brush-top jars.
- Advantages** Pro-Tite™ may be used on all plastic and metal pipe and fittings. It forms a resilient fluid and gas-tight seal, and is impervious to ammonia gas, aliphatic and aromatic petroleum distillates, including gasoline, chlorinated hydrocarbons, refrigerants, mineral and synthetic refrigeration oils, vegetable oil, fuel gasses, steam and water. Pro-Tite™ has a working temperature range of -100°F to 400°F. Pro-Tite™ sealed joints will withstand full vacuum and pressures up to 10,000 PSI.
- Directions** For best results always clean pipe threads before applying Pro-Tite™. Apply a thin uniform coating to the male threads of the joint, and assemble the joint. Tighten pipe fittings to manufacturer's recommended torque specifications. For best results allow 30 minutes for the sealant to set before returning the joint to service. Clean up with isopropyl or denatured alcohol solvent.
- Safety Data** Consult the *Material Safety Data Sheet* before using this product. Do not use this product until you have read all warning statements on the label and have consulted the MSDS. This product is only to be used by trained professional personnel. Use according to label directions. Wash thoroughly after using this product.
- Caution: Avoid eye and skin contact. Avoid breathing vapor. Do not take internally. Flammable, contains alcohol. Do not use or store near sources of ignition.
- First Aid: In case of eye contact, flush the affected area with water for 15 minutes. If irritation persists, get medical attention. For skin contact, wash thoroughly with soap and water. For ingestion, immediately drink 3-4 glasses of water. Do not induce vomiting. Get medical attention. For inhalation, remove the individual to fresh air if irritation occurs. If breathing does not return to normal within a few minutes, get medical attention.

Keep out of reach of children.



Package Metrics

- Producer of Goods** Specialty Chemical
- Domestic/Foreign Commodity** Domestic

SKU Specific Package Metrics

Catalog No.	Description	Size
CLL-1	Thread Sealant	1 1/3 oz. Squeeze Tube
CLL-4	Thread Sealant	4 oz. Brush Top Jar
CLL-16	Thread Sealant	16 oz. Brush Top Jar
UPC	Quantity per Case	
0095247086009	12	
0095247086016	12	
0095247086023	12	

Pro-Tite™ Thread Sealant